



X40T Thyroid Collar

- Standard size
- Large size
(over 20"/51cm neck)



X40H Half Aprons

- Velcro closure system
- 0.50 mm Pb equivalency

	<u>U.S. Size</u>	<u>Metric Size</u>
X40H-1	8x10"	20x25
X40H-2	12x14"	30x35
X40H-3	16x18"	40x45
X40H-4	18x24"	45x60
X40H-5	Set of 3 8x10" 12x14" 16x18"	Set of 3 20x25 30x35 40x45



X50H Humeral Shield

RADIATION PROTECTIVE EYEWEAR

All glasses are made with the highest quality SCHOTT glass and offer .75mm lead equivalency. Most styles are available in the following prescriptions: Single, Bifocal, and Progressive.



LT5050 Wrap

A lightweight fashionable wrap around nylon frame with a built in adjustable nose piece to distribute weight evenly over the nose. The ear-pieces are attached by spring hinges designed to withstand the rigors of hospital use.
*Weight: 85 grams



LT400 Fitover

The Fitover is designed to wear over prescription glasses. The oversized frame offers excellent protection and large viewing area.
*Weight: 92 grams



LT100 Ultralite

Manufactured of lightweight nylon and designed to provide protection for the entire eye area. The curved front and temple bar with rubber tips allow for a close secure fit. (Available in Black, Silver, Red, and Blue/Gray.)
*Weight: 58 grams



LT500 Maxi

The Maxi is slightly smaller than the Ultralite and includes a soft rubber nose bridge for increased comfort. (Available in Black, Silver, Red, and Navy)
*Weight: 63 grams



LT200 Euro

A standard safety frame which offers excellent protection and large viewing area. Universal form fit bridge.
*Weight: 75 grams



LT600 MetalFlex

The MetalFlex frames are constructed of optical steel base metal with either bronze or gold electroplate. Each pair is made with an adjustable non-slip silicone nose pad for maximum comfort.
*Weight: 62 grams



LT300 Astro

These widely used safety frames are quite flexible in temple length and the frame angle can be adjusted for a more custom fit. The front of the frame is fitted snugly to the brow to help alleviate fluid run-off into the eyes.
*Weight: 75 grams

***Weights on glasses are subject to change.