

# Protective Eyewear



# What makes Burlington Medical an outstanding medical protection partner?



Top quality products and more than 40 years of manufacturing experience.

---



Dedicated Account Managers who provide the highest level of customer service and product knowledge.

---



A vast assortment of comfortable, durable eyewear options for all your protective needs.

$\frac{20}{200}$

**E**

$\frac{200 \text{ FT}}{61 \text{ M}}$

$\frac{20}{100}$

**F P**

$\frac{100 \text{ FT}}{30.5 \text{ M}}$

$\frac{20}{70}$

**T O Z**

$\frac{70 \text{ FT}}{21.3 \text{ M}}$

$\frac{20}{50}$

**L P E D**

$\frac{50 \text{ FT}}{15.2 \text{ M}}$

$\frac{20}{40}$

**P E C F D**

$\frac{40 \text{ FT}}{12.2 \text{ M}}$

$\frac{20}{30}$

**E D F C Z P**

---

$\frac{30 \text{ FT}}{9.4 \text{ M}}$

$\frac{20}{25}$

**F E L O P Z D**

$\frac{25 \text{ FT}}{7.62 \text{ M}}$

$\frac{20}{20}$

**D E F P O T E C**

---

$\frac{20 \text{ FT}}{6.10 \text{ M}}$

$\frac{20}{15}$

**L E F O D P C T**

$\frac{15 \text{ FT}}{4.57 \text{ M}}$



**Your eyes are  
a critical tool  
for the work  
you do, and  
your vision is  
an important  
part of your life  
and well-being**

Scatter radiation and laser light both pose a serious health risk to your eyes, and every reasonable action should be taken to protect them. We offer an extensive line of eyewear that is designed for maximum comfort, weight distribution, and style. The variety of frames and features ensures there's a style to fit every need.

In this catalog, we present you with all the details you need to make informed decisions for the right eyewear. If you have any questions or would like to arrange an eyewear fitting, please reach out to your Account Manager or call +1-800-221-3466.

Model is wearing  
Stoggles in Lilac.  
See page 14 for  
additional color options.

## ANTI-FOG COATING

Our anti-fog coated lenses are essential for your radiation safety glasses. They offer an unparalleled level of fog resistance from a coating that will never fog when introduced into a humid environment. The anti-fog coating works on the molecular level to prevent vapors from cohering on the surface of the lens. Whether you are experiencing humidity from the room you're working in, or a bit of sweat, this coating will give you unclouded vision for the lifesaving tasks in your day. Only available with Non-prescription lenses or available in combination with anti-reflection coating on prescription lenses.

## ANTI-REFLECTIVE AND ANTI-FOG COATING

We've combined our anti-reflective and highly effective anti-fog coatings to create a highly transparent lens for procedures. This thin layer, applied to the lens surface, delivers both enhanced glare and fog prevention needed for procedures. This dual anti-reflective and anti-fog coating is available for both plano and prescriptions lenses.

## FRAME ENGRAVING

Personalize your glasses with the engraved frames to protect your investment. Engraving can be done on many of our frames so you can easily identify which glasses belong to you. There are limitations on the number of characters each frame can hold. If your request doesn't fit, we will work with you to find a solution.

## LENSES

We offer both plano (non-prescription) as well as prescription options for single vision, bifocal and progressive bifocal prescriptions. All front lenses come in 0.75mm Pb Lead protection.

## SIZING RECOMMENDATIONS

Each frame has detailed specifications of the frame size listed. Measurements are in millimeters (mm) unless otherwise stated.

There is no substitute for trying frames on, but it is helpful to narrow it down to the best options for your needs.

## SIDE SHIELDS

Side shields provide additional protection and coverage. Side shields are available for purchase to fit select frames. Frames that come with integrated side shields can be upgraded to 0.5mm protection. If leaded protection is not added, the side shields will be clear and offer no protection. Speak with an Account Manager for more details.

### KEY

Single Vision SV	Bifocal BIF
Anti-Fog AF	
Progressive PRG	Anti-Reflective & Anti-Fog AR/AF

All glasses come in a hard case with a retaining strap, cleaning cloth, and care instructions.

Please note that the weight of our eyewear products may vary within +/- 5% of the listed weight.

# OPTIONS ON EYEWEAR



SV	BIF
AF	
PRG	AR/ AF

**Ray-Ban RB4264 Ledged Glasses**  
**RAY-BAN RB4264 | WIDE FIT**

The classic Ray-Ban RB4264 is a square shape with a full rim frame made from lightweight plastic and adjustable nose pads. This pair of lenses comes with top-quality ledged glass lenses with a protection level of .75mm Pb lead.

● Black

Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
66g	57mm	18mm	145mm	46mm	132mm

Sphere: -7.00/+3.50    Cylinder: -5.50/+2.50  
Face Shape: Round, Square, Rectangle  
Side Shield Option: No



SV	BIF
AF	
PRG	AR/ AF

**Ray-Ban 7228 Ledged Glasses**  
**RAY-BAN 7228 | MEDIUM/AVERAGE FIT**

The Ray-Ban RB7228 is crafted with durable plastic construction designed for enduring wear and boasting the iconic Wayfarer shape great for a round face. This pair of lenses comes with top-quality ledged glass lenses with a protection level of .75mm Pb lead.

● Matte Black    ● Polished Opal  
Dark Grey

Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
75g	55mm	21mm	145mm	45mm	131mm

Sphere: -7.00/+5.00    Cylinder: -5.50/+5.00  
Detachable Amourx Side Protection Available: YES



SV	BIF
AF	
PRG	AR/ AF

**Ray-Ban Justin Ledged Glasses**  
**RAY-BAN JUSTIN | MEDIUM/AVERAGE FIT**

Ray-Ban Justin frames are Inspired by the iconic Wayfarer Sunglasses, this style has been redesigned with slightly larger rectangle lenses, and a rubberized frame finish to really give it that extra edge. This pair of lenses comes with top-quality ledged glass lenses with a protection level of .75mm Pb lead.

● Black

Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
68g	54mm	16mm	140mm	44mm	126mm

Sphere: -8.00/+5.00    Cylinder: -5.00/+5.00  
Side Shield Option: No



SV	BIF
AF	
PRG	AR/ AF


**Nike 5002 Ledged Glasses**  
**NIKE 5002 | SMALL/NARROW FIT**

This Lightweight injected plastic frame comes with cam action hinges snap that into place ensuring they stay on your head making them a great option for fast-paced environments. These lenses comes with top-quality ledged glass lenses with a protection level of .75mm Pb lead.

● Black    ● Gray

Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
53g	48mm	15mm	130mm	36mm	141mm

Sphere: -8.50/+4.00    Cylinder: -5.50/+5.00  
Detachable Amourx Side Protection Available: YES



## Nike Brazen Boost Ledged Glasses NIKE BRAZEN BOOST | MEDIUM/AVERAGE FIT


The lightweight plastic construction of these radiation protection eyewear, feature soft rubber nose pads and temple tips to provide extra comfort. This pair of lenses comes with top-quality ledged glass lenses with a protection level of .75mm Pb lead.



Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
63g	57mm	16mm	130mm	39mm	160mm

Sphere: -2.00/+2.00    Cylinder: -2.00/+2.00  
Face Shape: Round, Square, Oval  
Detachable Amourx Side Protection Available: YES

### KEY




All glasses come in a hard case with a retaining strap, cleaning cloth, and care instructions.

Please note that the weight of our eyewear products may vary within +/- 5% of the listed weight.




## Nike Brazen Ledged Glasses NIKE BRAZEN | WIDE FIT

A sleek design brought to you by Nike features a wraparound rectangular frame made of high-quality TR-90 Nylon. This pair of lenses comes with top-quality ledged glass lenses with a protection level of .75mm Pb lead.



Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
68g	60mm	16mm	135mm	37mm	136mm

Sphere: -2.00/+2.00    Cylinder: -2.00/+2.50  
Face Shape: Long, Round, Oval  
Side Shield Option: No



## Nike Brazen Fury Ledged Glasses NIKE BRAZEN FURY | WIDE FIT

This full-rim, rectangular shape lightweight plastic frame offers protection while providing premium levels of comfort. The rubber nose bridge ensures minimal slip making. Lenses comes with top-quality ledged glass lenses with a protection level of .75mm Pb lead.



Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
69g	60mm	17mm	130mm	43mm	172mm

Sphere: -2.00/+2.00    Cylinder: -2.00/+2.00  
Face Shape: Round, Square, Oval  
Detachable Amourx Side Protection Available: YES

# OPTIONS ON EYEWEAR




**Nike 7118 Ledged Glasses**  
NIKE 7118 | MEDIUM/AVERAGE FIT

Nike 7118 is a full rim frame for men and women that is made of plastic. This model features a square shape with a single bridge. This pair of lenses comes with top-quality ledged glass lenses with a protection level of .75mm Pb lead.



Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
57g	55mm	14mm	140mm	36mm	154mm

Sphere: -10.00/+4.00 Cylinder: -5.50/+5.50  
Face Shape: Long, Heart, Diamond  
Detachable Amourx Side Protection Available: YES




**Nike 7119 Ledged Glasses**  
NIKE 7119 | MEDIUM/AVERAGE FIT

Designed with premium injected material made from plant-based resin, which reduces impact on the environment. The inside temple tips feature flow lines for stability and grip. Lenses come with top-quality ledged glass lenses with a protection level of .75mm Pb lead.



Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
50g	53mm	17mm	140mm	31mm	146mm

Sphere: -7.00/+4.00 Cylinder: -5.50/+5.00  
Detachable Amourx Side Protection Available: YES




**Nike 7246 Ledged Glasses**  
NIKE 7246 | SMALL/AVERAGE FIT

The Nike 7246 features an acetate frame design for premium quality and durability. The frame comes with solid temple colorations with pop color accents and sneaker tread pattern. It includes Burlington's distortion-free SF-6 Schott glass lenses which offer a protection level of .75mm Pb lead equivalency.



Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
67g	54mm	17mm	140mm	38mm	155mm

Sphere: -10.00/+4.50 Cylinder: -5.50/+5.00  
Face Shape: Heart, Diamond  
Detachable Amourx Side Protection Available: YES




**Nike Premier Ledged Glasses**  
NIKE PREMIER | WIDE FIT

This frame offers a comfortable, lightweight frame with advanced optics, that provides precise vision from all angles and includes a rubber inlay nose pad and temple tips to provide extra comfort. They come with top-quality ledged glass lenses with a protection level of .75mm Pb lead.



Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
73g	60mm	13mm	135mm	39.3mm	170mm

Sphere: -3.00/+2.50 Cylinder: -3.00/+2.50  
Face Shape: Round, Square, Oval  
Detachable Amourx Side Protection Available: YES



## Nike City Icon Ledged Glasses NIKE CITY ICON | MEDIUM/AVERAGE FIT


The Nike City Icon Radiation Glasses feature rectangular frames made from Injected plastic material and made from at least 40% castor bean oil. This pair of lenses comes with top-quality ledged glass lenses with a protection level of .75mm Pb lead.



Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
66g	56mm	16mm	140mm	41mm	163mm


Sphere: -5.50/+3.50    Cylinder: -5.50/+3.00  
Side Shield Option: No

### KEY



All glasses come in a hard case with a retaining strap, cleaning cloth, and care instructions.

Please note that the weight of our eyewear products may vary within +/- 5% of the listed weight.



## Kate Spade Luella Ledged Glasses KS LUELLA | SMALL FIT

The Kate Spade Luella is a full rim that gives flexibility. The frame is acetate-made, thereby giving it numerous benefits such as light weight, durability, and secure fit. This pair of lenses comes with top-quality ledged glass lenses with a protection level of .75mm Pb lead.



Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
62g	51mm	18mm	140mm	38mm	143.3mm

Sphere: -10.00/+4.00    Cylinder: -5.50/+5.00  
Detachable Amourx Side Protection Available: YES



## Kate Spade Kattalin Ledged Glasses KS KATTALIN | MEDIUM/AVERAGE FIT

The Kate Spade Kattalin offers a rectangle shape design and feature a full rim acetate frame. This pair of lenses comes with top-quality ledged glass lenses with a protection level of .75mm Pb lead.



Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
61g	54mm	14mm	140mm	40mm	122mm

Sphere: -10.00/+4.00    Cylinder: -5.50/+5.00  
Detachable Amourx Side Protection Available: YES

# OPTIONS ON EYEWEAR



SV	BIF
AF	
PRG	AR/ AF

## Kate Spade Atalina Ledged Glasses KS ATALINA | MEDIUM/AVERAGE FIT

The rounded edges and upper corners give this frame a vintage feel. Crafted in lightweight, lustrous plastic, they're comfortable to wear. This pair of lenses comes with top-quality leaded glass lenses with a protection level of .75mm Pb lead.



Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
68g	49mm	16mm	140mm	39mm	145mm

Sphere: -10.00/+5.50    Cylinder: -5.50/+5.00  
Face Shape: Long, Heart, Diamond  
Detachable Amourx Side Protection Available: YES



SV	BIF
AF	
PRG	AR/ AF

## Under Armour Glacial Ledged Glasses UNDER ARMOUR GLACIAL | WIDE FIT

Under Armour's Glacial frames offer side shields for extra protection and non-slip pads to keep them on your face. This pair of lenses comes with top-quality leaded glass lenses with a protection level of .75mm Pb lead.



Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
78g	58mm	16mm	145mm	45mm	132mm

Sphere: -3.00/+3.00    Cylinder: -3.00/+3.00  
Face Shape: Oval, Rectangular  
Side Splash Shields Included, Option to Add 0.5mm LEV



AF	AR/ AF
----	-----------

## Under Armour Attack 2 Ledged Glasses UA ATTACK 2 | WIDE LARGE FIT

The Attack 2 features a durable injected plastic frame with a laser-etched trim, built-in rubber nose pads to prevent hair snagging, and TR90 hypoallergenic frames that meet ANSI Z87 standards. Auto-lock stop hinge ensures a secure fit. This pair comes with top-quality leaded glass lenses with a protection level of .75mm Pb lead.



Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
89g	63mm	16mm	130mm	45mm	142mm

Sphere: Non-RX Only    Cylinder: Non-RX Only  
Face Shape: Square, Rectangular  
Side Shield Option: No



SV	BIF
AF	
PRG	AR/ AF

## Under Armour 5043 Ledged Glasses UNDER ARMOUR 5043 | MEDIUM/AVERAGE FIT

These frames are made of plastic, which offers a combination of strength and lightweight comfort. The sustainable structure is made from Tritan Renew, a revolutionary recycled plastic that enhances their longevity. This pair of lenses comes with top-quality leaded glass lenses with a protection level of .75mm Pb lead.



Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
62g	55mm	15mm	140mm	40mm	125mm

Sphere: -6.00/+4.00    Cylinder: -5.50/+5.00  
Side Shield Option: No



AF | AR/  
AF

## Oakley Litebeam Ledged Glasses OAKLEY LITEBEAM | MEDIUM/AVERAGE FIT


They have a square shape with a durable O Matter, a nylon composite material that is designed to withstand everyday wear and tear. They also have adjustable wire-core temples and four interchangeable sizes of TruBridge nose pads. This pair of lenses comes with top-quality leaded glass lenses with a protection level of .75mm Pb lead.

 Satin Grey Smoke

Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
64g	55mm	18mm	140mm	40mm	160.5mm

Sphere: -4.00/+3.00    Cylinder: -4.00/+2.50  
Detachable Amourx Side Protection Available: YES

### KEY



All glasses come in a hard case with a retaining strap, cleaning cloth, and care instructions.

Please note that the weight of our eyewear products may vary within +/- 5% of the listed weight.



AF | AR/  
AF

## Oakley Holbrook Ledged Glasses OAKLEY HOLBROOK | MEDIUM/AVERAGE FIT

The Oakley Holbrook offers a timeless, classic design that is fused with modern Oakley technology in the O Matter, nylon composite material frame. This pair of lenses comes with top-quality leaded glass lenses with a protection level of .75mm Pb lead.

 Black     Steel

Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
79g	55mm	18mm	132mm	43mm	128mm

Sphere: Non-RX Only    Cylinder: Non-RX Only  
Face Shape: Round, Oval, Rectangular  
Side Shield Option: No



SV | BIF  
AF  
PRG | AR/  
AF

## Liberty Maxx 30 Ledged Glasses LIBERTY MAXX 30 | SMALL/MEDIUM FIT

These radiation protection glasses offer a stylish wrap-around frame constructed of high impact TR-90 nylon and feature rubberized soft grip temple bars and rubberized saddle bridge nose pad for extra comfort. This pair of lenses comes with top-quality leaded glass lenses with a protection level of .75mm Pb lead.

 Black

	Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
Reg:	73g	53mm	17mm	130mm	33mm	123mm
Lrg:	77g	55mm	22mm	130mm	39mm	133mm

Sphere: -10.00/+5.00    Cylinder: -5.50/+3.00  
Face Shape: Heart, Diamond  
Side Shield Option: No

# OPTIONS ON EYEWEAR



SV	BIF
AF	
PRG	AR/ AF

## Liberty X8-100S Ledged Glasses LIBERTY X8-100S | MEDIUM/AVERAGE FIT

The X8-100S has integrated side shields with the option to add protective shielding and is a purpose-built safety frame with ANSI Z87 rated protection. It's designed using Nylo-Flex™ material specifically built for ruggedness and flexibility. This pair comes with top-quality ledged glass lenses with a protection level of .75mm Pb lead.



Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
75g	54mm	16mm	130mm	35mm	124mm

Sphere: -7.00/+4.00    Cylinder: -5.00/+5.00  
Face Shape: Oval, Heart, Diamond  
Side Shield Option: Yes



SV	BIF
AF	
PRG	AR/ AF

## Pentax Classic 56 Ledged Glasses PENTAX CLASSIC 56 | MEDIUM/AVERAGE FIT

The Pentax Classic 56 Prescription Safety Glasses is a plastic frame that feature permanent brow guard side shields. This pair of lenses comes with top-quality ledged glass lenses with a protection level of .75mm Pb lead.



Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
73g	55mm	17mm	137mm	39.7mm	127mm

Sphere: -7.00/+4.00    Cylinder: -5.50/+5.00  
Face Shape: Square, Rectangular, Oval  
Clear Side Splash Shields Included, Option to Add 0.5mm LEV.  
Protection will be black and not clear.



SV	BIF
AF	
PRG	AR/ AF

## Bomber Clutch Ledged Glasses BOMBER CLUTCH | MEDIUM/AVERAGE FIT

These groundbreaking CLUTCH glasses combine modern aesthetics with tough safety features. Cam hinges and patented foam lining ensure comfort and security, while removable side shields provide ANSI Z87+ certified protection. This pair comes with top-quality .75mm Pb ledged glass lenses.



Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
81g	57mm	19mm	140mm	45mm	154mm

Sphere: -7.00/+5.00    Cylinder: -5.00/+5.00  
Face Shape: Square, Diamond  
Clear Side Splash Shields Included, Option to Add 0.5mm LEV.  
Protection will be black and not clear.



SV	BIF
AF	
PRG	AR/ AF

## Pentax Classic 10 Ledged Glasses PENTAX CLASSIC 10 | MEDIUM/AVERAGE FIT

The Pentax Classic 10 is a plastic frame with removable side shields that meet the ANSI Z87.1 industry standards. This pair of lenses comes with top-quality ledged glass lenses with a protection level of .75mm Pb lead.



Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
66g	52mm	16mm	140mm	38.7mm	120mm

Sphere: -10.00/+4.00    Cylinder: -5.50/+5.00  
Face Shape: Square, Rectangular  
Clear Side Splash Shields Included, Option to Add 0.5mm LEV.  
Protection will be black and not clear.



AF | AR/  
AF

## Jonathan Paul Cool Classic Ledged Glasses COOL CLASSIC | WIDE FIT | FITOVER FRAME


The Cool Classic is a flattering large fitover frame. Frames are lightweight Nylon with exclusive Scout Fit Temple design doesn't need adjustment and avoids added pressure on the ears. This pair of lenses comes with top-quality leaded glass lenses with a protection level of .75mm Pb lead.

● Black

Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
80g	59mm	17mm	140mm	45mm	135mm

Sphere: Non-RX Only    Cylinder: Non-RX Only  
Side Shield Option: No

### KEY



All glasses come in a hard case with a retaining strap, cleaning cloth, and care instructions.

Please note that the weight of our eyewear products may vary within +/- 5% of the listed weight.



AF | AR/  
AF

## Haven LSND Ledged Glasses HAVEN LSND | WIDE (FITOVER FRAME)

The Haven LSND glasses feature high quality nylon frames and bend-to-fit temples to allow for a comfortable fit. Plus, with their fits-over design, Haven LSND's can worn either alone or over prescription eyewear. This pair of lenses comes with top-quality leaded glass lenses with a protection level of .75mm Pb lead.

● Black

Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
98g	60mm	15mm	139mm	50mm	145mm

Sphere: Non-RX Only    Cylinder: Non-RX Only  
Face Shape: Square, Round  
Includes Side Protection 0.5mm LEV.  
Sides will be black and not clear.



AF | AR/  
AF

## Leader Archer Ledged Glasses LEADER ARCHER | LARGE FITOVER FRAME

The ARCHER is a large Fitover style frame to be worn alone with with your normal prescription glasses.. This unisex frame is made of nylon with a square shape lens and features a full rim construction. This pair of lenses comes with top-quality leaded glass lenses with a protection level of .75mm Pb lead.

● Black

Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
86g	65mm	13mm	126mm	52mm	146mm

Sphere: Non-RX Only    Cylinder: Non-RX Only  
Face Shape: Square, Rectangular  
Includes Side Protection 0.5mm LEV.  
Sides will be black and not clear.

# OPTIONS ON EYEWEAR



AF | AR/AF

## Bollé Spicy Ledged Glasses BOLLE SPICY | MEDIUM/AVERAGE FIT

A trendy, unisex model of safety glasses with built-in, impact resistant side shields. They are light and comfortable, and come with integrated side shields. This pair of lenses comes with top-quality leaded glass lenses with a protection level of .75mm Pb lead.

 Black/Blue

Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
69g	53mm	19mm	138mm	40mm	125mm

Sphere: -3.00/+2.50    Cylinder: -3.00/+2.00  
Clear Side Splash Shields Included, Option to Add 0.5mm LEV.  
Protection will be black and not clear.



SV | BIF  
AF  
PRG | AR/AF

## ArmouRx 5004 Ledged Glasses ARMOURX 5004 | WIDE FIT

The AMOURX 5004 combines a minimal design and staple color, making it an easy-to-wear safety frame. An integrated side shield and versatile shape provide a simple, sophisticated look suitable for all faces. This pair comes with top-quality leaded glass lenses with a protection level of .75mm Pb lead.

 Grey

Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
78g	60mm	13mm	135mm	40mm	133mm

Sphere: -3.00/+2.50    Cylinder: -3.00/+2.00  
Clear Side Splash Shields Included, Option to Add 0.5mm LEV.  
Protection will be black and not clear.



SV | BIF  
AF  
PRG | AR/AF

## Stoggles Rectangle Ledged Glasses STOGGLES RECTANGLE | SMALL/AVERAGE FIT

It's business meets barrier, keeping you sophisticated and safe at the same time. This polycarbonate frame has the option to add protection for the built in side shields. This pair of lenses comes with top-quality leaded glass lenses with a protection level of .75mm Pb lead.

 Denim     Moss     Lilac     Sangria

Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
64g	49mm	20mm	145mm	44mm	135mm

Sphere: -8.00/+4.00    Cylinder: -5.50/+3.00  
Face Shape: Rectangular  
Clear Side Splash Shields Included, Option to Add 0.5mm LEV.  
Protection will be black and not clear.



SV | BIF  
AF  
PRG | AR/AF

## ArmouRx 6019 Ledged Glasses ARMOURX 6019 | MEDIUM/AVERAGE FIT


ArmouRx 6019 glasses are lightweight, non-conductive, and meet ANSI Z87.1 standards. Made of high-quality plastic, they feature integrated side protection, rubberized nose pads, and adjustable rubberized temples for a superior fit. This pair comes with top-quality leaded glass lenses with a protection level of .75mm Pb lead.

 Gray/Orange

Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
80g	61mm	15mm	129mm	42mm	129mm

Sphere: -2.50/+2.50    Cylinder: -2.50/+2.50  
Face Shape: Square, Oval, Rectangular, Round  
Clear Side Splash Shields Included, Option to Add 0.5mm LEV.  
Protection will be black and not clear.

# OPTIONS ON EYEWEAR



## On Guard OG220S Ledged Glasses ON GUARD OG220S | WIDE/AVERAGE FIT

OnGuard 220S glasses are comfortable, lightweight, and made from non-conductive TR-Nylon. They feature a rubberized nose bridge for a superior fit and a dust dam for enhanced protection against dust and debris. This pair comes with top-quality ledged glass lenses with a protection level of .75mm Pb lead.



Black

Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
72g	57mm	15mm	135mm	25mm	155mm

Sphere: -5.00/+3.50    Cylinder: -5.00/+3.50  
 Face Shape: Round, Oval  
 Clear Side Splash Shields Included, Option to Add 0.5mm LEV.  
 Protection will be black and not clear.



## Pentax ZT500 Ledged Glasses PENTAX ZT500 | MEDIUM/WIDE FIT

The Pentax ZT500 are safety glasses that are ANSI-tested and have a durable front with short side shields, flexible temples, and an air seal around the eyes. They also have a removable head strap for a more secure fit. This pair of lenses comes with top-quality ledged glass lenses with a protection level of .75mm Pb lead.




Black/Blue

Weight (grams)	Lens Width	Bridge Width	Temple Length	Lens Height	Overall Width
77g	58mm	18mm	131mm	37mm	134mm

Sphere: -2.50/+2.50    Cylinder: -2.50/+2.50  
 Clear Side Splash Shields Included, Option to Add 0.5mm LEV.  
 Protection will be black and not clear.

### KEY



All glasses come in a hard case with a retaining strap, cleaning cloth, and care instructions.

Please note that the weight of our eyewear products may vary within +/- 5% of the listed weight.




## WARRANTY OF EYEWEAR:

1 year from the date of purchase against manufacturer defects (material and workmanship). This does not cover normal wear and tear, or negligent use of the product (including scratched lenses). Coating issues caused by cleaning lenses in a manner other than listed on the instructions for use will void the warranty. Repairs of eyewear are allowed for 2 years from the date of purchase. Repairs are limited to eyewear purchased from Burlington Medical. Prices and repair availability vary based on style and repairs being made.

Call your Account Manager today at 1-800-221-3466 or visit [BurMed.com](http://BurMed.com) for more details.

# Protection Eyewear at Your Fingertips

Shop our extensive line of eyewear from the comfort of your home, visit **BurMed.com**.



To place your order call your Account Manager at **1-800-221-3466**.



3 Elmhurst St, PO BOX 3194,  
Newport News, VA 23603